EE 535 Lab 2: Optical Measurements of Thin Films

Jonathan Hess

Abstract

The precise thickness of a thin absorbing wafer can be found using the transmission observed from a spectrophotometer. This procedure was done for amorphous silicon (aSi) and zinc selenide (ZnSe) wafers with unknown widths.

Introduction

Two wafers of different materials were observed with a spectrophotometer to yield measurements of The method evaluated in this paper to calculate the wafer film's width uses the transmition measured as well as the peaks and valleys of the resulting graph.

Experimental

During this lab measurements of two wafers were taken. These included reflection, transmission, absorbance, refractive index, as well as instrument error. One wafer was amorphous silicon and the other was of zinc selenide (ZnSe). What was not given or know was the thickness of these wafers. These measurements were done with the Cary 5000. This spectrophotometer has a range of 175-3300nm [5] and allows for the comparison of two samples at the same time. This was used to compare the sample wafer and a blank glass simultaneously. This allows to remove influence of the glass on the wafer's measurements.

The medium in the spectrophotometer is air so there will be an assumption that the refraction index will be exactly 1. There will also be an assumption that the wafer's thickness is uniform. This is because the graph in figure reflects a uniform thickness[3].

Definitions

Transmission (T) is the amount of light and electromagnetic radiation that passes through a media. $https://github.com/Jetsama/EE535/tree/main/Lab2 \\ Reflection~(R)~is~the~amount~of~light~and~electroscopic for the contraction of the contraction$ magnetic radiation that changes direction.

> Absorbance (A) is the measure of how much light is absorbed by a substance at a particular wavelength.

Theory

Interference-free transmission can be calculated with using following equations.

$$T_s = \frac{(1-R)^2}{1-R^2} \tag{1}$$

$$R = ((n_1 - n_2)/(n_1 + n_2))^2$$
 (2)

The equation 2 is that standard reflection coefficient for light striking a boundary between two different mediums. Where n_1 is the refraction index of first media and n_2 is the refraction index of the second media. Because the lab uses air as a medium this will be set to 1. Also because the equation is squared the order of these two refraction indices does matter. This derives equation 3 where s = the refraction index of the wafer.

$$R = ((s-1)/(s+1))^2 \tag{3}$$

Replacing R in equation 1 with equation 3 yields equation 4.

$$T_s = ((s-1)/(s+1))^2$$
 (4)

Solving for s gives refraction index in terms of transmission in equation 5.

$$s = \frac{1}{T_s} + (\frac{1}{T_s^2} - 1)^{1/2} \tag{5}$$

The equation for calculating interference fringe is

$$2nd = m\lambda \tag{6}$$

where m represents the integer for maximas and half integer for minima and λ .

$$T = \frac{Ax}{(B - Cx\cos(\phi)) + Dx^2} \tag{7}$$

$$A = 16n^2s \tag{8a}$$

$$B = (n+1)^3 (n+s^2)$$
 (8b)

$$C = 2(n^2 - 1)(n^2 - s^2)$$
 (8c)

$$D = (n-1)^3 (n-s^2)$$
 (8d)

$$\phi = \frac{\pi nd}{\lambda} \tag{8e}$$

$$x = exp(-\alpha d) \tag{8f}$$

Where $\lambda =$ wavelength, x = absorbance, d = width, s = refractive index, and n = complex refractive index

This equation is a simplification of the optical system equation where the extinction coefficient, k, has been set to 0. This is valid over most of the spectrum[9]. At the wavelengths tested silicon has a small k constant. This can be seem in figure where the value of the extinction coefficient rests at less than 0.05 during the measurement wavelengths.




Figure 1: Graph showing wavelength (nm) versus the extinction coefficient of silicon[10]

The local maxima and minima are defined by the 9 and 10 equations. Where T_M is the local maximums and T_m is the local minimums. These equations were derived from finding where $cos(\phi) = 1$ (where the T would be at a relative max) for the maximums and $cos(\phi) = -1$ for minimums. We can then replace $cos(\phi)$ from the with the +1 and -1 to

gain the T_M and T_m equations[9].

$$T_M = \frac{Ax}{(B - Cx) + Dx^2}\tag{9}$$

$$T_m = \frac{Ax}{(B+Cx)+Dx^2} \tag{10}$$

Using these maxima (or minima) the thickness of the thin film wafer can be found

$$d = \frac{\lambda_1 \lambda_2}{2(\lambda_1 n_{e2} - \lambda_2 n_{e1})} \tag{11}$$

To find the distance two extrema (either maximums or minimums) can be used along with equation 11. This was done for all extrema for all 3 wafers.

These equations equations were derived and then used on outputted data files were using a python script which can be found in the repository.

Results

After running the spectrometer there were data points for the two wafers. These data points included transmission, absorption and reflection index.

The amorphous silicon wafer transmission is shown in figure . This can be seen to be uniform thickness by the cosine wave imbued into the points. To use equation 11 extrema must be located. These were located and plotted on the graph as well as located in the table in figure .




Figure 2: Graph showing wavelength (nm) versus the transmission of Wafer aSi

The local minimums and maximums were calculated using the python module scipy. The values can be seen in figure below. Using these values combined with equation 8 the width of the wafer.

| num | maximums (nm) | minimums (nm) |
|-----|---------------|---------------|
| 1 | 659 | 667 |
| 2 | 691 | 705 |
| 3 | 731 | 751 |
| 4 | 781 | 807 |
| 5 | 841 | 875 |
| 6 | 921 | 965 |
| 7 | 1022 | 1077 |
| 8 | 1154 | 1228 |
| 9 | 1328 | 1435 |
| 10 | 1579 | 1730 |
| 11 | 1954 | 2213 |

Figure 3: Table with wavelength of maxima and minima for the aSi wafer

Using equation 11 a simple

The same process was done for the zinc




Figure 4: Graph showing wavelength (nm) versus the transmission of the ZnSe Wafer

Appendix

For information on the pure data or computational scripts there is a repository for this and other labs. https://github.com/Jetsama/EE535

References

- [1] Solid State Electronic Devices (2006, Prentice Hall), Ben Streetman, Sanjay Banerjee
- [2] An Entry Level Guide to Understanding Lasers (2008), Chapter 9.2.4, Jeff Hecht, 3rd ed.
- [3] R Swanepoel 1983 J. Phys. E: Sci. Instrum. 16 1214
- [4] Optical Processes in Semiconductors, Jacques I. Pankove
- [5] Cary 100/300/4000/5000/6000i/7000 Spectrophotometers User's Guide
- [6] How To Correctly Determine the Band Gap Energy of Modified Semiconductor Photocatalysts Based on UV-Vis Spectra, Patrycja Makuła, Michał Pacia, and Wojciech Macyk https://pubs.acs.org/doi/10.1021/acs.jpclett.8b02892
- [7] Silicon Photo Multipliers Detectors Operating in Geiger Regime: an Unlimited Device for Future Applications, Giancarlo Barbarino, Riccardo de Asmundis, Gianfranca De Rosa, Carlos Maximiliano Mollo, Stefano Russo and Daniele Vivolo
- [8] text
- [9] E Marquez et al 1992 J. Phys. D: Appl. Phys. 25 535
- [10] Handbook of Optical Constants of Solids, Edward D. Palik. Academic Press, Boston, 1985